



Forever Medical

## Portable Dental X Ray Machine YJ-DXP01



### Packing List


X-Ray Camera \*1

Charger \*1

User Manual \*1



Forever Medical



## Parameter

Type	Wireless Portable x-Ray Camra
Input Power-X ray unit	70 KVAC,5 mA
Input Power-Li-ion Charger	100-220VAC,50-60Hz,1.0A
Output Power-X ray unit	80w
Output Power-Li-ion battery capacity	14.8VDC,10A
Tube battery capacity	950mAh
Frequency	400KHZ
Exposure time	0.01-2.0S
Nominal focal spot value	0.4mm
Total filtration	1.5mmAl
Display mode	Digital
Weight	2.2kg
Packing Size	35x35x16cm

## Feature

High frequency inverter technology output is more stable  
Customs ultra-high capacity lithium battery is more durable  
One-click RVG digital model conversion more accurately  
Boot self-checking conventional fault automatic report errors  
Built-in temperature control,high voltage and low voltage protection safety

## PRODUCT SIZE



### Multiple usage scenarios

Can be held and mounted tripod



Note: This product does not include a tripod